



MWCCS Visit Code LDMS Reference Guide

Summary:

As of V105 visit (starting 01-Oct-2024), the MWCCS project employs a set of new Visit Codes within the Specimen Management Module of LDMS to denote different types of collections. When filling out the Visit Field labs are now required to implement the new codes in place of the previously used *Vst* code. The new Visit Codes are as follows:

LDMS VISIT UNIT CODE	DESCRIPTION
B/L	Baseline
CFU	Core Follow Up
SHR	Short
PFT	PFT
РК	РК
ОМВ	Oral MB
CRW	Crosswalk

The new visit codes can be implemented by 2 different ways:

- 1. Visit Specific Quick Add Templates: Each lab can create 7 unique Quick Add templates, one for each Visit Code. For example, a Baseline template can be created that includes only the specimens collected at the Baseline visit. When using this Baseline template the Visit Code will automatically be set to *B/L* and the set of specimens created in LDMS will be only those required for a Baseline visit.
- Single Quick Add Template + Manual Visit Code Selection: Each lab can establish only one Quick Add template with no default Visit Code that includes the specimens collected at every visit type. When completing data entry in LDMS the user must carefully select the appropriate Visit Code and manually remove specimens not collected at the applicable visit.





Method 1: Visit Specific Quick Add Templates

With the Visit Specific method, each lab must create 7 unique templates, one for each Visit Code within the *Quick Add Templates* utility of the Specimen Management Module in LDMS. When using these unique templates to add samples into the database, technicians will first choose the appropriate visit template within the *Quick Add* utility of LDMS. The template will automatically load the appropriate Visit Code and limited specimen collections expected for this visit on the Quick Add screen.

- 1. Navigate to *Quick Add* from the Specimen Management drop-down.
- 2. Apply the Project filter "MWCCS" to narrow down template list results.

I DMC	Specimen Management Storage Storage Reports Labels
	Test Specimen Management sion Administration Online Resources
	Quick Add
	Quick Add Templates Quick Add
Template Select a template 💌	Template Comments
Clear Template	Participant Information
Template Filters	Project 🗸 🗸
2. Mwccs	ID1 Add ID1 Upload ID1 File
Study	
Filter by ID2	OPIDs Add OPID
Visit Filter by visit value Visit Unit	OPID
Filter by visit unit	Enrollment Information
Clear Filters	ID2

3. Select the appropriate template from the Template drop down menu.



 Select Yes to confirm and load the template onto the Quick Add screen. The visit code and expected collections will auto-populate.

Confirm	х				
Would you like to load the template data? This will clear any data that you have previously entered.					
Yes No					





5. Fill-in or edit fields according to typical workflow.

Project		MWCCS		- *									
ID1 / PID				- * Add	d ID1 / PID	Upload	I ID1 / PID File						
OPIDs				Add C	OPID								
		OPID											Action
Enrollmont	Information												
Enroinnenu	mormation												
ID2 / Study	y	MWCCS		-									
Visit Informa	ation												
ID3				-									
Clinic		X826		-									
Collection I	Date	dd/MMM/y	'YYY	*									
Visit Value													
Visit Units		SHR		•	- Auto	o-pop	oulated vis	it code		Visit-sp	ecific co	ollectio	ns
									L				
 Primary Info 	ormation												
									*				Add New
#	Primary Type	Additive Type	Condition	Collection Time	Received Da	ate	Received Time	Volume	Volume Units	Additional Time	Additional Time Units	Other Spec	
1	BLD 👻	SST 🔻	SAT 🔻	HH:mm	dd/MMM/yyy	y –	HH:mm	8.5	ML	*	v		Edit 💌
Aligueto for													
- Aliquots for	Primary #1												
- Aliquots for	Primary #1—												Add New
Total Aliquots	Primary #1 Derivative Type	Sub A/D Type	Condition	Volume	Volume Units	Other S	Specimen ID						Add New





Method 2: Single Quick Add Template + Manual Visit Code Selection

With this method labs maintain only one Quick Add template that includes all collection types and has no default Visit Code. Technicians will first choose this template within the *Quick Add* utility of LDMS. The template will load a comprehensive set of specimen collections and will default Visit Code to null/empty. The technician must then carefully select the appropriate Visit Code from the visit code drop-down and manually remove any specimens not collected at the applicable visit.

- 1. Navigate to *Quick Add* from the Specimen Management drop-down.
- 2. Apply the Project filter "MWCCS" to narrow down template list results.
- 3. Select the appropriate template from the Template drop down. For example, *MWCCS V104 POS FEMALE*.

Template	
MWCCS V104	POS FEMALE 🗸
	Clear Template
Template Filt	ers
Project	
MWCCS	-
Study	
Filter by	ID2 🔻
Visit	
Filter by	visit value
Visit Unit	
Filter by	visit unit 🔻
	Clear Filters

4. Select *Yes* to confirm and load the template onto the Quick Add screen. The visit code will <u>not</u> auto-populate. The expected collections will auto-populate.

Confirm	х
Would you like to load the template data? This will clear any data that you have previously entered.	
Yes No	





5. Carefully input the appropriate visit code from the provided list of new codes in the *Visit Unit* field. For example, *SHR*, or Short Visit. Do <u>not</u> input the previously used *Vst* code.

-Participant Information																
Project		MWCCS			* *											
ID1 / PID		Add ID1 / PID Upload ID1 / PID File														
OPIDs		Add OPID														
		OPID											Acti	on		
- Enrollment Information-																
ID2 / Study		MWCCS			- *											
Visit Information																
ID3																
Clinic		X857														
Collection Date		dd/MMM/	<i>YYYY</i>		- *								_			
Visit Value						Careful	lv	input the	approp	riate v	isit	code.				
Visit Units		SHR			-	1	1	Dono	ot use. V	st.						
									<u></u>							
Primary Information																
															Add	New
# Primary Type	A Ty	dditive ype	Conditio	on	Collection Time	Received Date		Received Time	Volume	Volume Units		Additional Time	Additional Time Units	Other Spec		
1 BLD	- S	ST 👻	SAT	Ŧ	HH:mm	dd/MMM/yyyy	Ŧ	HH:mm	35	ML	Ŧ			201	Edit	-
2 URN	- N	ION -	SAT	-	HH : mm	dd/MMM/yyyy	Ŧ	HH:mm	15	ML	Ŧ		-	211	Edit	-
3 BLD	▼ E	DT 🔻	SAT	*	HH:mm	dd/MMM/yyyy	•	HH:mm	30	ML	-		-	203	Edit	-
4 VAG	- N	ION 👻	SAT	-	HH:mm	dd/MMM/yyyy	Ŧ	HH:mm	1	EA	•			212	Edit	-
5 BLD	▼ E	DT 👻	SAT	-	HH:mm	dd/MMM/yyyy	Ŧ	HH:mm	2	CRD	Ŧ		-	224	Edit	-

6. Fill-in or edit all remaining fields according to typical workflow including removing all specimen rows not applicable to this visit type.