

## Storage Setup Examples

**LDMS User Support** 

Phone: 716-834-0900 x7311 E-mail: ldmshelp@fstrf.org www.frontierscience.org/ldms

**LDMS Training Team** 

Phone: 716-834-0900 x7466

E-mail:

ldms.training@fstrf.org

## Chest Freezer

The chest freezer holds racks in a linear configuration. In the example below each rack is uniform but multilple capacities may be used. Racks may be numbered (or named) on top to easily identify when the lid is opened. These may be applied in the Main View tab after adding the template.


Label: Chest  Type: Freezer   ▼	Freezer							
Temperature: -70 C  Configuration Design:  4x4 rack 1  4x4 rack 2  4x4 rack 3  4x4 rack 4  4x4 rack 5  4x4 rack 5  4x4 rack 6  4x4 rack 7  4x4 rack 8	1	2	3	4	5	6	7	8

## LN 2 Tank no divisions

In the example provided, the 1x20 racks have been added to the LDMS in a linear setup. In reality the racks are in a tank in a staggered configuration due to the circular nature of the container (or on a carousel). The number is placed on top of the rack to be easily identified from above.

Label: Nitrogen Freezer  Type: Freezer  Temperature: Liquid Nitrogen  ▼	
Configuration Design: 1x20 rack 1	_
1x20 rack 3 1x20 rack 4 1x20 rack 5 1x20 rack 6 1x20 rack 7 1x20 rack 8 1x20 rack 9 1x20 rack 10 1x20 rack 11 1x20 rack 12	
1x20 rack 13 1x20 rack 14 1x20 rack 15 1x20 rack 16 1x20 rack 17 1x20 rack 18 1x20 rack 20	

## LN2 Tank with divisions



In this example the circular divisions need to be flattened. This is done by creating a starting point. Label the device if not yet labeled. The 9 racks in the diagram are in a staggered configuration. In the LDMS each sector should be formatted to a square (ie. 3x3)

Freezer	Sector A	Sector B	Sector C	Sector D	Sector E	Sector F							
		Soctor A											
	Sector A												
Shelf	Rack 1	Rack 2	Rack 3		Box 1	1 x8 Rack							
	Rack 4	Rack 5	Rack 6	1	Box 2								
	Rack 7	Rack 8	Rack 9	1	Box3								
				-	Box 4								
					Box5								
					Box 6								
					Box 7								
					Box8								