


Quick Reference Guide for Specimen Management in LDMS for the Web

LDMS User Support

Phone: 716-834-0900 x7311 E-mail: ldmshelp@fstrf.org www.fstrf.org/ldms LDMS User Manual

LDMS Training Team Phone: 716-834-0900 x7466 E-mail: ldms.training@fstrf.org

Filters are available on several pages and are used to narrow down information that is displayed, to locate desired records.



Specimen Records |< Participant 5 of 8 🔻 REPRIEVE Project Edit Participant 💌 0012345L PID A specimen record in LDMS for the Web records OPIDS specimens collected at each visit for a specific participant PROTOCOL A5332 Edit PROTOCOL -Visits for PROTOCOL A5332 Collection Date SID Clinic Action Visit 21/Apr/2016 NOSID Edit 🔻 Primary Specimens for Visit 0 Ent. 21/Apr/2016 Additive Type Available Volume Other Specimen Global Specimen ID Status Collection Time Primary Type Action SST SAT 5 ML 9005-001MMC00-000 07:00 BLD Edit 🔻 9005 -001MNG00-000 07:00 BLD SST SAT 20 MI Edit 🔻 SAT 20 ML 9005-001MPD00-000 💮 🔚 🖂 🦺 07:00 BLD EDT Edit 🔻 5 ML 9005-001MQK00-000 💮 🔚 🖂 📣 SPO SAT 07:00 BLD Edit 🔻 9005-001MRG00-000 💮 🚟 🖂 🚜 07:00 URN NON SAT 1 ML Edit 🔻 07:00 SAT 5 ML 9005-001MSF00-000 💮 🚟 🖂 🦺 BLD SST Edit -9005-001MTC00-000 💮 🔚 🖂 🧀 20 ML 07:00 BLD SST SAT Edit 🕶 9005-001MVB00-000 💮 🚟 🖂 🦚 20 ML Edit 🔻 Edit -Aliquots for 9005-001MMC00-000 Global Specimen ID Status Available Volume Other Specimen Specimen Action SAT 9005-001MMC00-001 N/A 1.8 ML SER Edit 🕶

The Fundamentals

Specimen Processing Information should be entered for all REPRIEVE samples.

Click on the **Edit** button for each Primary and Aliquot. Be sure to enter the following:

Primary


- Processing Date
- Processing Time
- •Processing Tech Initials Aliquot
- •Frozen Date
- •Frozen Time

LDMS Condition Codes should be used to indicate when expected specimens are not collected or specimen issues are encountered during processing. Common condition codes* include the following: SNC (Sample Not Collected), QNS (Quantity Not Sufficient) and SHV (Short Volume).

Condition code updates should be entered at the primary or aliquot level, in the **Sample Condition** field.

In addition, the *Reasons Sample Not Collected* field in the Primary and Aliquot Edit windows of the LDMS can be utilized to document why expected samples are not available.

*A full listing of condition codes is available in the Reports module, Miscellaneous category or in the user manual. Specimen Labels can be printed via the **Edit** button at the Visit level.



be used, if printing on previously used sheets of labels.

Entering Specimens: Quick Add Templates

Template A5332 ENTRY - NON-NETWORK SITES FSTRF-defined Quick Add Templates available for REPRIEVE protocols make it easier for laboratory staff to enter samples into LDMS for the Web, by pre-Clear Template populating information into the Quick Add screen of Specimen Management. When the desired template is selected, specimen information is populated based Template Filters on the protocol-defined collections and required processing for study visits. The Project user must then complete the required information for record completion: **ID1**, REPRIEVE Clinic, ID3, Draw Date, Draw Time, Received Date, and Received Time, ID2 and make any needed adjustments to volumes and condition codes, prior to Visit adding the record. Visit Unit Quick Add Clear Filters Template Comments Ship to BRI per assigned schedule. Document Processing Start and Freeze times Participant Information Project REPRIEVE PID 0111111C OPIDs Add OPID OPID Action Enrollment Information PROTOCOL A5332 * Visit Information SID NOSID Clinic 101 06/May/2016 Collection Date Visit Value 0 Visit Units Ent Primary Information Add New **Additive Type** Volume **Volume Units** Primary Type **Collection Time Received Date** Received Time BLD SST 05:00 06/May/2016 06:00 ML 5 Edit 🔻 2 BLD SST 05:00 06/May/2016 06:00 20 ML Edit 3 BLD EDT 05:00 06/May/2016 06:00 20 ML Edit 4 BLD SP0 05:00 06/May/2016 06:00 5 ML Edit 5 URN NON 05:00 06/May/2016 06:00 1 ML Edit 🔻 > Aliquots for Primary #1 Add New **Total Aliquots Derivative Type** Sub A/D Type Volume **Volume Units** Other Specimen ID SER N/A ▼ 1.8 ML Edit ▼