

LDMS Specimen Management Codes and Volume for MTN 001

Use the table below as a guide when entering 001 specimens in the LDMS. For each test, enter the corresponding codes and volume (primary and aliquot) in the LDMS Specimen Management module.

Test	Primary	Additive	Derivative	Sub Add/Der	Primary Volume	Aliquot Volume	Units
Vaginal Biopsies	VGL	NON	TIS	N/A	Variable	1	N/A
PBMC for Intracellular Tenofovir	BLD	CPS	CEL	MET	Variable	1.0	ml
Cervicovaginal Lavage	CVL	NSL	CVL	N/A	Variable	1.0	ml
Plasma for Storage	BLD	EDT	PL1/2	N/A	Variable	1.0–2.0	ml
Serum for Tenofovir	BLD	NON	SER	N/A	Variable	1.0–2.0	ml
Cytology Brushes	CER	PBS	CTB	N/A	Variable	1.0	N/A