

## LDMS Specimen Management Codes and Volume for MTN 002

Use the table below as a guide when entering 002 specimens in the LDMS. For each test, enter the corresponding codes and volume (primary and aliquot) in the LDMS Specimen Management module.

Test	Primary	Additive	Derivative	Sub Add/Der	Primary Volume	Aliquot Volume	Units
Blood for Tenofovir Levels	BLD	NON	SER	N/A	4.0	1.5	ml
Amniotic Fluid	AMN	NON	AMN	N/A	2.0–4.0	1.0–2.0	ml
Cord Blood	CRD	NON	SER	N/A	4.0	1.5	ml
Placental tissue	PLC	NON	PLC	N/A	1.0–2.0	1.0–2.0	cm
Endometrial tissue	END	NON	END	N/A	1.0–2.0	1.0–2.0	cm