

LDMS Specimen Management Codes and Volume for MTN 003B

Use the table below as a guide when entering 003B specimens in the LDMS. For each test, enter the corresponding codes and volume (primary and aliquot) in the LDMS Specimen Management module.

Test	Primary	Primary Additive	Primary Volume	Primary Units	Aliquot Derivative	Aliquot Sub Add/Der	Aliquot Volume	Aliquot Units
Serum for Storage	BLD	NON	15.0	ml	SER	N/A	2.0	ml
Urine for Storage	URN	NON	15.0–60.0	ml	URN	N/A	Variable	ml