

Patient, Administration & Labels Quiz

1. True or False: If you wish to change patient information for all records logged into your LDMS under a defined patient ID, group, and protocol/study, then this change should be made in the Patient Module.

2. The Backup Tracking module allows you to track backups performed at your laboratory. Which of the following statements are true?
 - A. You can add, modify and delete backup entries
 - B. The Set Range button allows you to change how your previously-entered backup entries are displayed.
 - C. You can view all entries, all backups from the past X days, or all backups that have been performed since a particular date.
 - D. The Backup Tracking module in the LDMS does not back up your data. This feature is solely for the purpose of tracking backups that have already been performed at your lab.

3. To generate a label in the Labels module, the following pieces of information must be defined prior to clicking the Execute button:
 - A. Group
 - B. Group and Label Type
 - C. Group, Label Type and all of the of the following: Received Date, Received Batch Number, LDMS Specimen Number or Global Specimen ID
 - D. Group, Label Type and at least one of the following: Received Date, Received Batch Number, LDMS Specimen Number or Global Specimen ID

4. True or False: The User Configuration module, within the Administration tab, can be used to create or assign an individual user name and password for each staff member at your lab?

5. Your LDMS automatically times out after it has sat idle for 10 minutes, and you are then forced to log back into the system. You wish to change the idle timeout to 30 minutes. Which module in the Administration tab should be accessed to change the idle timeout?
 - A. System Configuration
 - B. User Configuration
 - C. Instrument Setup
 - D. Group Configuration

6. True or False: The LDMS forces you to use the same printer for generating both reports and labels.

7. As laboratory manager, you wish to restrict your laboratory staff to allow view (read-only) permission in the Specimen Management and Storage modules, while allowing for full permissions (view, add, delete and modify) in the Shipping module. Which module in the Administration tab would you use to setup these types of restrictions?

- A. System Configuration
- B. User Configuration
- C. Instrument Setup
- D. Group Configuration

8. True or False: You can use the Group Configuration module in the Administration tab to make changes to pre-defined government group identifiers (government groups include ACTG, WIHS, CHAVI, HPTN, VTN).

9. True or False: The LDMS supports various type of labels. It is possible to select and set a default label size and format within your local LDMS.

10. True or False: The LDMS allows a user to skip spaces within a sheet of labels and begin printing labels from a particular position and forward.

11. True or False: It is possible for a user to generate (print) custom label formats for government groups.
