



LDMS User Support
Phone: 716-834-0900 x7311
E-mail: ldmshelp@fstrf.org
www.ldms.org/contact/

LDMS Training Team
Phone: 716-834-0900 x7466
E-mail: ldms.training@fstrf.org

The Reagent Logging module is used to track reagents used in the laboratory and uses a defined set of reagent codes. The module will record the Lot no., expiration date, manufacturer and start/end dates. Data entered into the reagent log is exported to Frontier Science.

Reagent Logging

	Reagent Name	Reagent Lot #	Lot Exp. Date	Manufacturer	Start Date	End Date
1	DMS	123456	11/Oct/2017	FISHER SCIENTIFIC	01/Mar/2017	01/Aug/2017
2	DMS	321654	17/Jan/2018	FISHER SCIENTIFIC	01/Aug/2017	09/Nov/2017
3	DMS	445566	09/May/2018	FISHER SCIENTIFIC	09/Nov/2017	
4						

1	Reagent Name	Reagent Lot #	Lot Exp. Date	Manufacturer	Start Date	End Date
	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Go to **Administration > Reagent Logging**

1. Below the grid complete the data fields. The following fields are required:

Reagent Name, Reagent Lot #, Lot Exp. Date, and Manufacturer

Enter the **Start Date** (optional) and **End Date** (if available)

2. Add the entry to the log by clicking the **+** button

An entry can be amended or corrected by highlighting the entry in the log. Add or change the information in the fields and save the changes by clicking the **Save** button 