

LDMS Reference Guide

Shipping Module: QA/QC Report

The Shipping QA/QC Report lists the QA/QC status of the specimens in the shipment. This report is intended to be printed and included with the other manifest reports in the physical shipment to the receiving lab.

Pending Shipments

Shipment Number	Shipment Format	Destination Lab	Shipment Temperature	Setup Date	Shipment Date	QA/QC	
16	LDMS			19/Aug/2021		Not Performed	Edit/Ship

1 - 1 of 1 results

- Generate Manifest
- Shipment Storage Report
- Shipping Container Report
- Assign Tests
- Print Labels
- QA/QC
- QA/QC report
- Delete

Shipping QA/QC Report

Box #1-1

Global Spec ID	Scanned Global Spec ID	Status	QA/QC Tech Initials	QA Date
9999-0002SC00-001	Manual	Passed	AZ	19/Aug/2021
9999-0002SC00-002	Manual	Failed	AZ	19/Aug/2021
9999-0002SC00-003	9999-0002SC00-003	Passed	AZ	19/Aug/2021
9999-0002SC00-004	9999-0002SC00-003	Failed	AZ	19/Aug/2021
9999-0002SC00-005		not performed		
9999-0002SC00-006		not performed		

Manual QA/QC

Barcode Scanner QA/QC

Select a container: Box #1-1

Container 1 of 1

OK	!	OK	!	?	?	?	?	?
?								